# Your Port Access **Advantage\***

---

*Implanted Port Education Program*

---

03 INTRODUCTION

---

03 WHAT IS A BARD\* PORT

---

03 ABOUT DUAL LUMEN PORTS

---

04 HOW IS YOUR PORT PLACED

---

04 HOW IS YOUR PORT USED

---

05 HOW IS YOUR PORT ACCESSED

---

05 CARE OF YOUR PORT

---

06 WHAT IS A POWERPORT\* DEVICE

---

06 CECT SCANS / POWER INJECTION

---

06 IDENTIFICATION

---

07 IMPORTANT SAFETY INFORMATION

---

08 TELL YOUR CLINICIAN

---

08 COMMON QUESTIONS AND ANSWERS

---

09 NOTES FOR QUESTIONS TO ASK

## INTRODUCTION

Your doctor has recommended a Bard Access Systems port because your treatment requires frequent delivery of medications into your bloodstream. With your port, you'll be able to receive IV treatment without repeated needle sticks in your arm or leg veins.

Please read all of the information in this Patient Guide. It is designed to help you understand and feel comfortable with your port.

This booklet is a reference. Your doctor and nurse are still the best source of information. You should talk to them about any questions or concerns that you may have. Always follow their directions for the use and care of your port. You can also visit [www.portadvantage.com](http://www.portadvantage.com) for more information about ports from Bard Access Systems.

## WHAT IS A BARD\* PORT?

Your Bard port is a small device placed beneath your skin. It is a small basin that is sealed with a rubber top called a septum. The port connects to a small, soft tube called a catheter. The catheter is placed inside one of the large central veins that take blood to your heart. When a special needle is put into the septum, it creates "access" to your bloodstream. Medicine and fluids can be given through the needle and blood samples withdrawn.

Ports come in a variety of shapes and sizes. Your doctor has chosen the port that is best for your body size and therapy needs.

## ABOUT DUAL LUMEN PORTS

Various treatments may be needed along with chemotherapy. For example, chemo patients may need blood therapies, antibiotics, anti nausea drugs, pain medicines and hydration. A dual lumen port may be ideal for such patients. Also a dual lumen port may be indicated for power injection of contrast media and offers the benefits of providing for power-injected CT scans.

## HOW IS YOUR PORT PLACED?

The port is placed during a short procedure. The procedure is most often done under local anesthetics. You may be able to go home the same day. Your doctor will place the port just under the skin through a small cut. He will then place the end of the catheter into your blood vessels. The port will be placed in the best position for your treatment. Often, the port is placed on the upper chest just below the collar bone. You will see a small raised area beneath the skin where the port is placed.

## HOW IS YOUR PORT USED?

Your doctor or nurse can give medicines or fluids or take blood samples through your port. They will not need to stick your arm veins with a needle. This makes it more comfortable for you. Because the port places medicines into the large central veins, the medicines mix better in the blood. The medicines are diluted so they are less harmful to your veins.

The port can be used to deliver medicine in two ways. One way, called “bolus injection”, delivers medicine all at once with the needle in place for only a short time. Another way is a slow delivery called a “continuous infusion”. For this type of medicine, a dressing will be placed over the needle. The dressing will hold the needle in place for a longer period of time. A small tube will connect the needle to the medicine. The port also provides an easy way to collect blood samples. When a blood sample is needed a special needle is placed into the port. Blood is withdrawn through the needle.

## HOW IS YOUR PORT ACCESSED?

Your doctor or nurse will use the port when they need to give medicine or fluids or withdraw blood samples. To do this, they will access the port by placing a special needle into the port. You may feel a pricking during needle insertion. This sensation often gets milder over time. Ask your doctor or nurse what you should expect.

A=Safety Infusion Set


B=Septum

C=Port

D=Skin

E=Catheter

F=Blood Vessel

## CARE OF YOUR PORT.

After your port is placed, your doctor will place a small bandage over the wound. For the first few days, you should avoid heavy exertion and follow any special guidance from your doctor to care for the small incision. Once the incision has healed, it will not take any special care and you can resume normal activities.

### **About heparin locks**

Sometimes blood can clot and block the catheter. This clot stops medicines and fluid from flowing through it. Blood will not clot when it is thinned with a medicine called heparin. To prevent clots, ports with open-ended catheters are often filled with a mixture of saline and heparin. This process is called a heparin lock. If your port will not be used for a long period of time, the doctor or nurse will usually change the heparin lock every four weeks.

If you are allergic to heparin or have Heparin Induced Thrombocytopenia (HIT), you should inform your doctor or nurse of this anytime the port is used.

### **About saline locks**

Another option to a heparin lock is the “saline lock.” Ports with a valved catheter, called a Groshong\* Catheter, are often filled with saline. Saline locks simplify port maintenance. Also saline locks help patients avoid allergic reactions to heparin. If your port will not be used for a long period of time, the doctor or nurse will usually change the saline lock every four weeks.

## WHAT IS A POWERPORT\* DEVICE?

Power injectable ports provide access for power-injected Contrast Enhanced Computer Tomography (CECT) or CT scans. CT scans produce quick, accurate images of your body. CT scans using a power injectable port will help the medical team better manage your care without having to stick your arm with a needle.

## CECT SCANS / POWER INJECTION

Contrast Enhanced Computer Tomography (CECT) scans are simple, safe and non-invasive procedures that aid in disease diagnosis and treatment. They are more sensitive than x-rays and allow the doctor to see small differences in your soft tissue that they may not see with normal x-rays.

Before doing a CECT scan, the clinician will inject a special fluid through your port. This fluid makes the picture taken with the scan clearer. For best results, a fluid called contrast media is delivered into your blood at a high rate. This process is called power injection.

A PowerPort\* device when used with a special needle allows clinicians to perform power-injected CECT scans without putting a needle into your arm or leg veins.

## IDENTIFICATION

Some special features help identify a PowerPort\* device. The port has a special triangle shape and three bumps on top of the septum. The bumps are called palpation bumps. Your doctor or nurse can feel these bumps and the triangle port shape to know that you have a PowerPort\* device.

You should have received a special packet of information when a PowerPort\* device is placed. This packet includes an ID card, bracelet, and key ring that help identify your PowerPort\* device. Tell your clinician if you have a PowerPort\* device whenever your port is accessed for a procedure.

If you have a PowerPort\* device, a trained medical professional can identify the unique shape and design that distinguishes it from other ports. Tell your clinician that you have a PowerPort\* device before your port is accessed for a procedure.

A trained clinician can also recognize the PowerPort\* device in an X-ray, due to a special symbol located on the port, only visible under x-ray.

## IMPORTANT SAFETY INFORMATION

A port is not for everyone—especially patients with a history of forming blood clots, who have had previous vascular access surgery, or who are not emotionally prepared to have an implanted medical device. Like any vascular access procedure, there is always a risk of complications, including venous blood clots, skin erosion, infection, a collapsed lung, or clotting of the port catheter. Talk to your physician or nurse about these risks, and whether a port or other treatments are right for you.

### ***Who should not receive a port?***

Patients who:

- Have or suspected of having an infection
- Have a history of forming blood clots
- Have a body size that will not allow for proper port placement or port access
- Have had the port insertion site exposed to radiation
- Are not emotionally prepared to have an implanted medical device

### ***What are the possible complications associated with an implantable port?***

A port is not for everyone. It is important the you talk with your physician or nurse about the risks associated with a port and port use. Potential complications include:

- Blood clot formation
- Redness and inflammation at the port insertion site
- Skin erosion
- Bleeding at the insertion site
- Scarring at the insertion site
- Port catheter breakage
- Irregular heart beat
- Vein erosion
- Allergic reaction to the materials of the port
- Port rotation or flipping while implanted
- Leakage of medicine, IV fluid, or chemotherapy from the port
- Improper port access leading to medicine, IV fluid, or chemotherapy being injected into surrounding tissue of the port
- Port clotting

Use of an implanted port carries risks associated with a minor surgical procedure and vascular access. Potential complications include:

- Internal Bleeding
- Nerve damage
- Collapsed lung
- Fluid buildup around the lungs
- Blood clot formation
- Accidental cutting or puncturing of blood vessels

## TELL YOUR CLINICIAN.

As a new patient with a port, the best way to be involved is share information and concerns with your doctor or nurse. Speak up. Ask about your concerns, or if you see anything that seems odd.

What to report to your Doctor or Nurse

- You have an implanted port.
- If you have a Bard\* PowerPort\* device.
- If you notice any redness or swelling at the site of your port after the wound heals.
- If you have a fever.
- If you have allergies to medicines or materials.
- If you have an allergy to heparin.
- If you have heparin induced thrombocytopenia (HIT).
- If you have ever been given blood thinning medicines such as heparin or warfarin.
- If you have been treated with radiation before.
- If you have ever been diagnosed with lung disease.
- If you have ever been diagnosed with, or treated for blood clots in your veins (venous thrombosis).
- If you have ever been diagnosed with any tissue disease or suffered from tissue erosion.

- If you have been diagnosed with or tested for “pinch-off” syndrome.
- If other clinicians have ever had difficulty withdrawing blood or infusing fluids through your port. Let them know if other clinicians have asked you to change position to allow blood or fluid to flow.

## COMMON QUESTIONS AND ANSWERS

### ***Will the port affect my daily life?***

After the port has been placed and the small incision has healed, you should be able to return to your daily life. Ask your doctor or nurse about specific activities and when to resume them.

### ***Will I need to wear a bandage over the port?***

A bandage will be needed until the incision heals. Afterwards, a bandage is not required when the port is not in use. If you are receiving a treatment a bandage may be applied to secure the needle.

### ***Do I have to stop wearing certain types of clothes?***

Ask your doctor or nurse. The answer will depend on where your port is placed.

**Who pays for the port?**

Insurance policies vary. Check with your insurance company.

**Will the port set off security alarms?**

Security systems may detect the small amount of metal in the port. If this happens show your Patient Identification Card.

**How long will I have my port?**

The port can stay in place as long as your doctor decides that you need it.

**Can the device be removed if I no longer need it?**

Yes. The port can be removed in a simple procedure similar to the one used to place it.

**Can I get a CT procedure with a PowerPort\* device?**

Yes. The materials used in the PowerPort\* device are safe for use in CT and CECT procedures.

**Can I get an MRI procedure with a Port in place?**

Yes. You may get an MRI with this port. There may be some conditions, so be sure your clinician checks the MRI guidelines for the port.

**What do I do if I have a question or suspect a problem?**

After the procedure you will see some redness and tenderness in the area. This should go away in 24 to 48 hours. However, contact your doctor or nurse if you notice unusual changes in the skin such as increased swelling, redness or soreness. If you have pain, fever, chills, shortness of breath, or dizziness, contact your doctor right away.

**What do I do if my doctor or nurse is unavailable?**


You will be given a wallet-sized Medical Alert Identification Card. This card has important information about the port. With this information other medical personnel will be able to help you.

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**Bard Access Systems, Inc.**