

CyAn ADP

FL#	Excitation	Em Filter	Example Fluorophores
FL1	488	530/40	FITC, GFP, Alexa488, Rhodamim123, Newport Green, JC-1 (green), DiO
FL2	488	575/25	PE, *Cy3, *TRITC, TMRE, Auto Orange, JC-1 (red), *DiL, *DsRed
FL3	488	613/20	PI, PE-TxRed
FL4	488	680/30	PI, PE-Cy5, PerCP, 7AAD (For PerCP-Cy5.5 or PE-Cy5.5 use 700/20)
FL5	488	750LP	PE-Cy7
FL6	405	450/50	DAPI, Hoechst, Pacific Blue, CFP, BV421
FL7	405	530/40	Aqua Viability Dye (Pacific Orange use 580/30)
FL8	635	665/20	APC, Alexa647 (For Alexa700 use 700/20), Cy5
FL9	635	750LP	APC-Cy7, APC-e780
* Not optimal with a 488 excitation			